

Method 508 MOD.Date Extracted 2-20-02Page 96Ad#6CMS-17627

Lab #	Description	Sample vol/wt	Drip rate ml/min	Surr add	Spk add	IS add	Final vol.	T.V. pos.	pH	Count		
L.B.		500 ml		1 ml			1 ml	10B	7			
Ref 1					1 ml			10A				
Ref 2					1 ml			9A				
3825								9B				
3826	See 3/27							8A				
3827								7B				
3828								6A				
3829	F.B.							8B				
3830								7A				
3831								8A				
3832	F.B.							7B				
3823								7A				
3824								8B				
Remarks <u>Surrogate manually calculated</u>												
drip rate must be between 5-15 ml/min for liq-liq extractors												
Extraction type - <u>SEP. RUN.</u> Solv. lot# <u>V 42 EM</u> NaCl lot# <u>V39618</u> Na2SO4 lot# <u>V 46623</u>												
STANDARD	INT STD LOG#											
CALIBRATION	<u>HP05062</u>											
REF SPIKE	<u>2193</u>											
MDL/MDRF												
Q. CHECK:												
INT. STD.	<u>022162</u>											
SURROGATE	<u>021202</u>											
LPC												
ENDRIN/DDT												

Analyst A.D.
revision 4/12/01Verified by hyReviewed by DN
Date 3/29/02

Sequence Name: C:\HPCHEM\1\SEQUENCE\MAR2202.S

* Comment:

Operator:

Data Path: C:\HPCHEM\1\DATA\mar2202\

Pre-Seq Cmd:

Post-Seq Cmd:

Method Sections To Run

On A Barcode Mismatch

(X) Full Method

(X) Inject Anyway

() Reprocessing Only

() Don't Inject

Line Type	Vial	DataFile	Method	Sample Name
-----------	------	----------	--------	-------------

1 Sample	1	DFTPP	GCMS0308	
2 Sample	2	CAL40	GCMS0308	
3 Sample	29	5900	GCMS0308	gc/ms scan 3/13
4 Sample	30	5901	GCMS0308	gc/ms scan 3/13
5 Sample	3	LB	GCMS0308	gc/ms scan 2/20
6 Sample	4	REF1	GCMS0308	gc/ms scan 2/20
7 Sample	5	REF2	GCMS0308	gc/ms scan 2/20
8 Sample	6	3825	GCMS0308	gc/ms scan 2/20
9 Sample	7	3826	GCMS0308	gc/ms scan 2/20
10 Sample	8	3827	GCMS0308	gc/ms scan 2/20
11 Sample	9	3828	GCMS0308	gc/ms scan 2/20
12 Sample	10	3829	GCMS0308	gc/ms scan 2/20 fb
13 Sample	11	3830	GCMS0308	gc/ms scan 2/20
14 Sample	12	3831	GCMS0308	gc/ms scan 2/20
15 Sample	13	3832	GCMS0308	gc/ms scan 2/20 fb
16 Sample	14	3823	GCMS0308	gc/ms scan 2/20
17 Sample	15	3824	GCMS0308	gc/ms scan 2/20
18 Sample	1	DFTPP2	GCMS0308	
19 Sample	2	CAL40A	GCMS0308	
20 Sample	16	LBA	GCMS0308	gc/ms scan 2/25
21 Sample	17	REFA1	GCMS0308	gc/ms scan 2/25
22 Sample	18	REFA2	GCMS0308	gc/ms scan 2/25
23 Sample	19	4038	GCMS0308	gc/ms scan 2/25
24 Sample	20	4039	GCMS0308	gc/ms scan 2/25
25 Sample	21	4125	GCMS0308	gc/ms scan 2/25
26 Sample	22	4126	GCMS0308	gc/ms scan 2/25
27 Sample	23	4127	GCMS0308	gc/ms scan 2/25
28 Sample	24	4128	GCMS0308	gc/ms scan 2/25
29 Sample	25	4129	GCMS0308	gc/ms scan 2/25
30 Sample	26	4130	GCMS0308	gc/ms scan 2/25
31 Sample	27	4131	GCMS0308	gc/ms scan 2/25
32 Sample	28	4132	GCMS0308	gc/ms scan 2/25
33 Sample	31	LBB	GCMS0308	gc/ms scan 2/25
34 Sample	32	REFB	GCMS0308	gc/ms scan 2/25
35 Sample	33	REFB1	GCMS0308	gc/ms scan 2/25
36 Sample	34	4028	GCMS0308	gc/ms scan 2/25
37 Sample	35	4029	GCMS0308	gc/ms scan 2/25
38 Sample	36	4030	GCMS0308	gc/ms scan 2/25
39 Sample	37	4031	GCMS0308	gc/ms scan 2/25 fb
40 Sample	38	4032	GCMS0308	gc/ms scan 2/25
41 Sample	39	4033	GCMS0308	gc/ms scan 2/25
42 Sample	40	4034	GCMS0308	gc/ms scan 2/25


Last Modified: Mon Mar 25 08:04:26 2002

Page: 1

DFTPP

Data File : C:\HPCHEM\1\DATA\MAR2202\DFTPP.D
 Acq On : 22 Mar 2002 4:09 pm
 Sample :
 Misc :
 MS Integration Params: GOOFY.P
 Method : C:\HPCHEM\1\METHODS\GCMS0308.M (RTE Integrator)
 Title : 209 625/TCLP calibration table

Vial: 1
 Operator:
 Inst : GC/MS 209
 Multiplr: 1.00



Spectrum Information: Average of 26.836 to 26.873 min.

Target Mass	Rel. to Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Result Pass/Fail
51	198	30	60	35.0	78845	PASS
68	69	0.00	2	0.0	0	PASS
69	198	0.00	100	46.4	104595	PASS
70	69	0.00	2	0.0	0	PASS
127	198	40	60	43.1	97224	PASS
197	198	0.00	1	0.0	0	PASS
198	198	100	100	100.0	225445	PASS
199	198	5	9	6.6	14951	PASS
275	198	10	30	23.1	52146	PASS
365	198	1	100	2.4	5477	PASS
441	443	0.01	100	98.8	31464	PASS
442	198	40	110	73.8	166485	PASS
443	442	17	23	19.1	31848	PASS